# SALTER

## Maxview Analyser

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


---

Downloaded from [www.vandenborre.be](http://www.vandenborre.be)

## HOW DOES THE MAXVIEW ANALYSER SCALE WORK?

The Salter Maxview Analyser Scale uses BIA (Bio Impedance Analysis) technology which passes a tiny electrical impulse through the body to determine fat from lean tissue, the electrical impulse cannot be felt and is perfectly safe. Contact with the body is made via stainless steel pads on the platform of the scale.

This method simultaneously calculates your personal weight, body fat, total body water, BMI, Basal Metabolic Rate and muscle mass, giving you a more accurate reading of your overall health and fitness. This scale will also display weight, body fat and total body water for children and young adults aged 6-17 years. Muscle mass, BMI and BMR are not shown for this age range.

This scale stores the personal data of up to 12 users. As well as being an analyser scale, this scale can be used as a conventional scale.

## BODY FAT - WHAT DOES IT MEAN?

The human body is made up of, amongst other things, a percentage of fat. This is vital for a healthy, functioning body - it cushions joints and protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat can be damaging to your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror.

This is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measure of fitness than weight alone – the composition of your weight loss could mean you are losing muscle mass rather than fat - you could still have a high percentage of fat even when a scale indicates 'normal weight'.

## BODY WATER - WHY MEASURE IT?

Body water is the single most important component of body weight. It represents over half of your total weight and almost two thirds of your lean body mass (predominantly muscle).

Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body's temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be 'normal' but your body hydration level could be insufficient for healthy living.

## BASAL METABOLIC RATE (BMR) - WHAT DOES IT MEAN?

Your Basal Metabolic Rate (BMR) is an estimation of the energy (measured in calories) expended by the body, at rest, to maintain normal body functions. This continual work makes up about 60-70% of the calories the body uses, and includes the beating of your heart, respiration, and the maintenance of body temperature. Your BMR is influenced by a number of factors, including age, weight, height, gender, dieting and exercise habits.

## WHY SHOULD I KNOW MY MUSCLE MASS?

According to the American College of Sports Medicine (ACSM), lean muscle mass may decrease by nearly 50 percent between the ages of 20 and 90. If you don't do anything to replace that loss you're losing muscle and increasing fat. It is also important to know your muscle mass % during weight reduction. At rest, the body burns approximately 110 additional calories for each kilo of muscle gained. Some advantages of gaining muscle mass include:

- Reversing the decline in strength, bone density and muscle mass that accompanies age
- Maintenance of flexible joints
- Guide weight reduction when combined with a healthy diet.

**NORMAL HEALTHY RANGE OF BODY FAT & TOTAL BODY WATER PERCENTAGE**

MALE			FEMALE				
AGE	BODY FAT	BODY WATER	AGE	BODY FAT	BODY WATER		
6	12-20%	60%+	6	14-23%	50%+		
7	13-20%						
8	13-21%						
9	13-22%						
10	13-23%						
11	13-23%						
12	12-23%						
13	12-22%						
14	11-21%						
15	10-21%						
16	10-20%						
17	10-20%						
18-29	8-19%		60%	18-29		20-30%	50%
30-39	11-21%						
40-49	13-23%		55%	40-49		24-34%	47%
50-59	15-25%						
60-69	17-27%		50%	60-69		28-38%	45%

Figure 1

	BELOW NORMAL BODY WATER %	NORMAL & ABOVE BODY WATER %
LOW BODY FAT %	Fitness program recommended, ensure sufficient water intake	Fitness program recommended
HEALTHY BODY FAT %	Weight and body fat OK but insufficient water intake	Ideal balance body fat, muscle mass and weight
HIGH BODY FAT %	Change of diet required, ensure sufficient water intake	Change of diet required

Figure 2

For those using Athlete Mode: note that athletes may have a lower body fat range and a higher body water range than shown above, depending on their type of sport or activity.

Body water measurement results are influenced by the proportion of body fat and muscle. If the proportion of body fat is high, or the proportion of muscle is low then the body water results will tend to be low.

It is important to remember that measurements such as body weight, body fat and body water are tools for you to use as part of your healthy lifestyle. As short term fluctuations are normal, we suggest you chart your progress over time, rather than focus on a single days reading.


This product is for domestic use only. Pregnant women should only use the weight function.

**Professional medical guidance should always be sought before embarking on diet and exercise programmes.**

**CAUTION: Do not use this product if you have a pacemaker or other medical device fitted. If in doubt, consult your doctor.**

### WHAT DOES MY BMI VALUE MEAN?

The Salter Maxview Analyser Scale calculates your BMI value. Below is a useful chart to show the correlation between height, weight and the BMI categories Underweight, Normal, Overweight and Obese.

Your BMI category can be identified using the table below.

BMI Category	BMI Value Meaning	Health Risk Based Solely On BMI
< 18.5	Underweight	Moderate
18.6-24.9	Normal	Low
25-29.9	Overweight	Moderate
30+	Obese	High

**IMPORTANT:** If the chart/table shows that you are outside of a 'normal' healthy weight range, consult your doctor before taking any action.

BMI classifications are for adults (age 20+ years) only.

## QUESTIONS AND ANSWERS

### **How exactly is my body fat and water being measured?**

The Salter Maxview Analyser Scale uses a measurement method known as Bioelectrical Impedance Analysis (BIA). A minute current is sent through your body, via your feet and legs. This current flows easily through the lean muscular tissue, which has a high fluid content, but not easily through fat. Therefore, by measuring your body's impedance (i.e. its resistance to the current), the quantity of muscle can be determined. From this, the quantity of fat and water can then be estimated.

### **What is the value of the current passing through me when the measurement is taken? Is it safe?**

The current is less than 1mA, which is tiny and perfectly safe. You will not be able to feel it. Please note however, that this device should not be used by anyone with an internal electronic medical device, such as a pacemaker, as a precaution against disruption to that device.

### **If I measure my body fat and water at different times during the day, it can vary quite considerably. Which value is correct?**

Your body fat percentage reading varies with body water content, and the latter changes throughout the course of the day. There is no right or wrong time of the day to take a reading, but aim to take measurements at a regular time when you consider your body to be normally hydrated. Avoid taking readings after having a bath or sauna, following intensive exercise, or within 1-2 hours of drinking in quantity or eating a meal.

### **My friend has a Body Fat Analyser made by another manufacturer. When I used it I found that I got a different body fat reading to that on my Salter Maxview Analyser Scale. Why is this?**

Different Body Fat Analysers take measurements around different parts of the body and use different mathematic algorithms to calculate the percentage of body fat. The best advice is to not make comparisons from one device to another, but to use the same device each time to monitor any change.

### **How do I interpret my body fat and water percentage readings?**

Please refer to the Body Fat and Water tables included with the product. They will guide you as to whether your body fat and water reading falls into a healthy category (relative to your age and sex).

### **What should I do if my body fat reading is very 'high'?**

A sensible diet, fluid intake and exercise program can reduce your body fat percentage. Professional medical guidance should always be sought before embarking on such a program.

### **Why are the body fat percentage ranges for men and women so different?**

Women naturally carry a higher percentage of fat than men, because the make-up of the body is different being geared towards pregnancy, breastfeeding etc.

### **What should I do if my water percentage reading is 'low'?**

Ensure that you are regularly taking sufficient water and work towards moving your fat percentage into the healthy range.

### **Why should I avoid using the Salter Maxview Analyser Scale whilst pregnant?**

During pregnancy a woman's body composition changes considerably in order to support the developing child. Under these circumstances, body fat percentage readings could be inaccurate and misleading. Pregnant women should therefore only use the weight function.

## GETTING STARTED

### Step 1: Batteries

Insert the 4 x AAA batteries supplied into the battery compartment, ensuring the + and - terminals are correctly orientated.

### Step 2: Carpet Feet


For use on carpet please attach the carpet feet provided.

### Step 3: Set Time


- 1) Check scale is off.
- 2) Press and hold **SET** button until hour flashes on display. Adjust hour using **▲** and **▼** buttons (note clock is 24 hour). Press **SET** button to confirm.
- 3) Minutes flash. Adjust minutes using **▲** and **▼** buttons. Press **SET** button to confirm.

## PERSONAL DATA ENTRY


- 1) Press the **SET** button to turn the scale on.
- 2) While the user number is flashing, select a user number by pressing the **▲** or **▼** buttons. Press the **SET** button to confirm your selection.


- 3) The weight mode display will flash. Press **▲** button until the correct symbol is flashing, then press the **SET** button.

- 4) The height display will flash. Press **▲** or **▼** as necessary to set your height, then press the **SET** button.

- 5) The age display will flash.  
Press ▲ and ▼ as necessary to set your age, then press the **SET** button.


#### Child Mode:

When the age is set from 6-12 years child mode is automatically selected. A small child icon will appear, as illustrated below.

**Note:** In child mode only weight, body fat and body water are displayed.

- 6) The male or female symbol will flash.  
Press ▲ button until the correct symbol is flashing, then press the **SET** button.

#### Athlete Mode:

When the age is set from 15-50 years Athlete Mode is available.

To select Athlete Mode press the ▲ button until  is displayed, then press the **SET** button.

An athlete is defined as a person who is involved in intense physical activity of approximately 12 hours per week and who has a resting heart rate of approximately 60 beats per minute or less.

- 7) The display will show your settings, then switch off. The memory is set.  
8) Repeat procedure for a second user, or to change user details.


**Note:** To update or overwrite the memorised data, follow the same procedure, making changes as required.

**WEIGHT & BODY FAT READINGS**

Position the scale on a flat level surface.

- 1) Press the **SET** button.
- 2) Select your user number by pressing the **▲** or **▼** buttons.
- 3) **WAIT** while the display reconfirms your personal data then shows a zero reading.
- 4) When zero is displayed, **WITH BARE FEET**, step onto the platform and stand still. Ensure your feet are in good contact with the 4 metal pads on the platform. Your feet will cover the display.
- 5) After a few seconds you will hear a beep. Step off the scale.
- 6) Your weight will be displayed followed by your body fat percentage, total body water percentage, muscle mass percentage, BMI and Basal Metabolic Rate (BMR).

## WEIGHT READING ONLY

This scale features our convenient step-on operation. Once initialised the scale can be operated by simply stepping straight on the platform – no more waiting!

### Set Weight Mode

- 1) Press the platform centre and remove your foot. The time is displayed.
- 2) Press the platform centre and remove your foot. Zero is displayed.
- 3) Press **SET** button until kg, lb or st is displayed.
- 4) Wait for scale to switch off.

### Initialising Your Scale

- 1) Press the platform centre and remove your foot. The time is displayed.
- 2) Press the platform centre and remove your foot.
- 3) Wait for scale to switch off.

**This initialisation process must be repeated if the scale is moved. At all other times step straight on the scale.**

### Operating Your Scale

- 1) Step on and stand with your feet over the display.
- 2) When you hear a beep, step off.
- 3) Your weight is displayed.
- 4) Zero is displayed, then the time.
- 5) Scale switches off.

Tip: To display the time when the scale is off press the platform centre and remove your foot.

## WARNING INDICATORS

'Lo' = replace battery.

'O-Ld' = Weight exceeds maximum capacity.

'Err' = out of range, incorrect operation or poor feet contact.

## TECHNICAL SPECIFICATION

80 cm - 220 cm  
2'7" - 7'3"  
6 - 100 years  
Male/Female  
12 user memory  
Child mode = 6-12 years

Athlete mode = 15-50 years  
d = 0.1 % body fat  
d = 0.1 % body water  
d = 0.1 % muscle mass  
d = 1 kcal BMR

## ADVICE FOR USE AND CARE

- Always weigh yourself on the same scale placed on the same floor surface. Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- Placing your scale on a hard, even floor will ensure the greatest accuracy and consistency.
- It is recommended that you measure at the same time of the day, preferably early evening before a meal for the most consistent results.
- Your scale rounds up or down to the nearest increment. If you weigh yourself twice and get two different readings, your weight lies between the two.
- Clean your scale with a damp cloth. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics.
- Treat your scale with care - it is a precision instrument. Do not drop it or jump on it.

## TROUBLESHOOTING GUIDE

- You must have bare feet when using body fat and total body water measurements.
- The condition of your skin on the bottom of your feet can affect the reading. To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on the scale.

## GUARANTEE

This product is intended for domestic use only. Salter will repair or replace the product, or any part of this product, (excluding batteries) free of charge if within 15 years of the date of purchase, it can be shown to have failed through defective workmanship or materials. This guarantee covers working parts that affect the function of the scale. It does not cover cosmetic deterioration caused by fair wear and tear or damage caused by accident or misuse. Opening or taking apart the scale or its components will void the guarantee. Claims under guarantee must be supported by proof of purchase and be returned carriage paid to Salter (or local Salter appointed agent if outside the UK). Care should be taken in packing the scale so that it is not damaged while in transit. This undertaking is in addition to a consumer's statutory rights and does not affect those rights in any way. For UK Sales and Service contact HoMedics Group Ltd, PO Box 460, Tonbridge, Kent, TN9 9EW, UK. Helpline Tel No: (01732) 360783. Outside the UK contact your local Salter appointed agent.

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