Tensio

USER MANUAL | GUIDE D'UTILISATION

NL Handleiding | **IT** Manuale di istruzioni | **ES** Manual de instrucciones

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Terrailon®

www.terraillon.com

Downloaded from www.vandenborre.be

Made for / Compatibilité

-  iPhone® 4S/5/5C/5S/6/6Plus
-  iPod® Touch 5th generation
-  iPad® 3/4/Air
-  iPad® Mini
-  Android 4.3
-  Bluetooth Smart 4.0 / Bluetooth Smart Ready

iPhone is registered trademark of Apple Inc.

Android is a trademark of Google Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

Other trademarks and trade names are those of their respective owners.

Tensio

Downloaded from www.vandenborre.be

EN	TERRAILLON CONNECTED BLOOD PRESSURE MONITOR INSTRUCTION MANUAL	4
FR	TENSIOMÈTRE CONNECTÉE TERRAILLON NOTICE D'UTILISATION	8
NL	TERRAILLON POLSBLOEDDRUKMETER HANDLEIDING	12
IT	SFIGMOMANOMETRO TERRAILLON MANUALE DI ISTRUZIONI	16
ES	TENSIÓMETRO DE MUÑECA TERRAILLON MANUAL DE INSTRUCCIONES	20
DE	HANDGELENK-BLUTDRUCKMESSGERÄT TERRAILLON BEDIENUNGSANLEITUNG	24
PT	MONITOR DE PRESSÃO SANGÜINEA DE PULSO TERRAILLON MANUAL DE INSTRUÇÕES	28

Thank you very much for selecting the Terraillon Wrist Blood Pressure Monitor TENSIO.




This unit allows you to monitor your blood pressure. This device is not intended to be a diagnostic device. Always contact your doctor for advice, diagnostic and treatment.

FIND MORE INFORMATION ON

<http://www.terraillon.com>

PRODUCT OVERVIEW

- SYS** Systolic Blood Pressure
- DIA** Pression artérielle diastolique
- Pul min** Pulse
- ERROR** Error
- mmHg** Unit
- Bluetooth** Data transmitting
- Battery** Low Battery
- Memory** Mémoire
- Time** Time [Hour/Minute - Month/Day/Year]
- IHB** Irregular Heartbeat Detector

FIRST USE


- A** Press and hold "SET" button to set time. Press "MEM" button to change the numeral then press "SET" button to confirm. After confirming [TIME] and [DATE], the LCD will display "Done".
- B** Download the Terraillon Wellness Coach application.

- C** Create your account.


- D** Turn on Bluetooth on your Smartphone (Setting > Bluetooth > ON).
- E** Press and hold "START" to start pair-up.

START MEASURING

- A** Tie the cuff.
- B** Press "START" to turn on the monitor, and it will complete the measurement process automatically.

RECORDS

Press "MEM" button to access the memory. The monitor will display the calculated average of the last three readings first. Press "MEM/UP" button "SET/DOWN" button to rotate the history records. "MEM/UP" to go forward "SET/DOWN" to go backward.

TIPS FOR MEASUREMENT

	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	≥180
DIA	<80	80-84	85-89	90-99	100-109	≥110

* The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH).

TROUBLESHOOTING

- E2** The cuff is very tight.
- E3** The pressure of the cuff is excess.
- E9** Product has not been activated.
- E10/11** The monitor detected motion while measuring.
- E20** Measurement process does not detect the pulse signal.
- E21** Measure incorrectly.
- Eexx** A calibration error occurred.
- Lo** Low Battery. Recharge the product with USB cable.

WARRANTY & ENVIRONMENT PROTECTION

This product is guaranteed for 2 years against material and manufacturing defects. During this period, such defects will be repaired free of charge (proof of purchase must be presented if the scale is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence. If you have a complaint please contact the store where you purchased the product first.

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

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Safety Information

The below signs might be in the user manual, labeling of other components.

They are the requirement of standard and using.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
C € 0123	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION - Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "MANUFACTURER"		
SN	Symbol for "SERIAL NUMBER"		Symbol for "DIRECT CURRENT"
	The Bluetooth Combination Mark		Symbol for "Authorised Representative in the European Community"
	Symbol for "MANUFACTURE DATE"		Symbol for "Including RF transmitter"

CAUTION

- * It is intended for adult indoor use only. Pregnant women, pre-eclamptic patients and patients with severe obesity don't use the device. If you need, please consult professional doctor.
- * This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the wrist or for functions other than obtaining a blood pressure measurement.
- * Please use the device under specified environment by user manual, otherwise the accuracy of the device will be influenced.
- * Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Please start or end medical treatment basing solely on physician's treatment advice.
- * If you are taking medication, consult your physician to determine the most appropriate time for your measurement. Never change a prescribed medication without your physician's consent.
- * This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's wrist and fingers will become anaesthetic, swollen, and even purple due to a lack of blood.
- * If the pressure of the cuff exceeds 40 kPa (300 mmHg), the unit will automatically deflate. If the cuff don't deflate when its pressure exceeds 40 kPa (300 mmHg), detach the cuff from the wrist and press the START/ STOP button to stop inflation.
- * Do not use the monitor under the conditions of strong electromagnetic field (e.g. medical RF equipment) that radiates interference signal or electrical fast transient/ burst signal.
- * The maximum temperature that the applied part can be achieved is 42.5°C while the environmental temperature is 40°C.
- * The device is not AP/APG equipment. It is not suitable for use in the presence of a flammable anesthetic mixture with air (or oxygen, nitrous oxide).
- * Please keep the unit out of reach of infants, children or pets, since inhalation or swallowing of small parts is dangerous or even fatal.
- * Please use ACCESSORIES and detachable parts specified / authorised by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user / patient.
- * The patient is an intended operator. The patient can measure, transmit data and charge battery under normal circumstances and maintain the device and its accessories according to the user manual.
- * The blood pressure monitor, and the cuff are suitable for use within the patient environment. If you are allergic to dacron or plastic, please don't use this device.
- * The device is not intended for PATIENT transport outside a healthcare facility.
- * This device cannot be used with HF surgical equipment at the same time.
- * There is a PTC current limiter in the monitor, which specification is 8V and 0.5A. When the voltage and current exceed the limiting value, the monitor will stop working.
- * If Luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.
- * The device is not suitable for public use.
- * Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- * When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: too frequent and consecutive multiple measurements; the application of the CUFF and its pressurization on any wrist where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the wrist on the side of a mastectomy.
- * Do not apply the cuff over a wound, otherwise it can cause further injury.
- * Do not inflate the cuff on the same limb which other monitoring ME EQUIPMENT is applied around simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME EQUIPMENT.
- * Please make sure the unit functions safely and it is in proper working conditions before use. Don't service or maintain while the device is in use.
- * If you have any problems with this device, such as setting up, maintaining or using, please contact with SERVICE PERSONNEL of Terrailon. Don't open or repair the device by yourself.
- * Using it in case to result in prolonged impairment of the circulation of the blood of the PATIENT.
- * Please report to Terrailon if any unexpected operation or events occur.
- * Cleaning: Dust environment may affect the performance of the unit. Please use the soft cloth to remove the dirt of the device and cuff before and after use.
- * Calibration: The manufacturer does not require such preventive inspections or calibration by other persons and will make available on request of circuit diagrams, component part list, etc.
- * Disposal: Degraded sensors may result in inaccurate measurement while loosened electrodes may cause the monitor's failure to power on. Please dispose of ACCESSORIES, detachable parts, and ME EQUIPMENT according to local guidelines.

Tips for User

It can cause inaccuracy if the measurement is taken in the following circumstances.

1. Within 1 hour after dinner or drinking.
2. Immediate measurement after tea, coffee, smoking.
3. Within 20 minutes after taking a bath.
4. When talking or moving your fingers.
5. In a very cold environment.
6. When you want to discharge urine.

To obtain the best performance, please follow below instructions.

1. Put in a dry place and avoid the sunshine.
2. Avoid immersing it in the water. Clean it with a dry cloth in case.
3. Avoid shaking and collision.
4. Avoid dusty environment and unstable temperature surrounding.
5. Use the slightly damp cloth to remove the dirt.
6. Avoid washing the cuff.

CAUTION

1. The battery of Tensio is built-in rechargeable lithium-ion battery, please do not disassemble it by the unauthorized maintenance personel.
2. Under the normal using, it can charge power about 300 times, if the battery cannot charge the power normally or the blood pressure monitor cannot use normally, please connect with the authorized maintenance personel. If measured three times per day, and the battery is fully charged, it can be used for about 20 days.
3. Storge and use the blood pressure monitor at the cool, dry and ventilated environment. Avoid to approach to the fire and the heat source, or it will cause the battery explode.
4. Only can use the Teraillon's authorized USB cable to charge the power. You cannot use the blood pressure monitor during the process of charging.
5. During the process of charging, the blood pressure monitor display.  When the charging is finished, please pull the plug in time.
6. When charging, shall not touch charging connector and the patient simultaneously.

Positioning the Cuff

1. Remove all accessories (watch, bracelet, etc) from your left wrist. If your physician has diagnosed you with poor circulation in your left wrist, use your right wrist.
2. Roll or push up your sleeve to expose the skin.
3. Apply the cuff to your left wrist with your palm facing up.
4. Position the edge of the cuff about 1-2 cm.
5. Fasten the wrist cuff around your wrist, leaving no extra room between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.
6. Patients with Hypertension:

The middle of the cuff should be at the level of the right atrium of the heart; Before starting measurement, please sit comfortably with legs uncrossed, feet flat on the floor, back and arm supported.

- Resting for 5 minutes before measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same wrist, or as directed by a physician.

Irregular Heartbeat Detector

This Blood Pressure Monitor is equipped with an intelligent function of Irregular Heartbeat (IHB) Detector. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this equipment will light up the IHB symbol on the screen when displaying the measuring result.


CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies every in one day, it also affected by the way you tie your cuff and your measurement position, so please take the measurement at the same condition.
2. The varies of the pressure is greater if the person take medicine.
3. Waiting at least 3 minutes for another measurement.

Why the blood pressure I get from the hospital is different from home?

The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

The attention need to pay when you measure your blood pressure at home:

If the cuff is tied properly.

If the cuff is too hight or too loose.

If the cuff is tied on the wrist.

If you feel anxious pressured.

You had better take deep breath 2-3 times before beginning.

Advice: adjust yourself for 4-5 minutes until you calm down.

If the result is the same if measuring on the right wrist?

It is ok for both wrists, but there will be some different results for different person, so suggest you measure the same wrist every time.

Specification

Power supply	3.7V 420mAH Built-in rechargeable lithium ion battery, 5V / 1A USB
Display mode	Digital LCD V.A.46.5x36.5mm
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0kpa-40kpa (0mmHg-300mmHg) Measurement pressure: 5.3kPa-30.7kPa (40mmHg-230mmHg) pulse value : (40-199)beat/minute
Accuracy	Pressure: 5°C-40°C with ± 0.4 kpa (3mmHg) pulse value: $\pm 5\%$
Normal working condition	Temperature: 5°C to 40°C Relative humidity $\leq 85\%$ Atmospheric pressure: 86kPa to 106kPa
Storage & transportation condition	Temperature: -20°C to 60°C RH: 10% to 93% Atmospheric pressure: 50kPa to 106kPa
Measurement perimeter of the wrist	About 13.5cm-21.5cm
Weight	Approx.110g
External dimensions	Approx. 79.8×72.5×13.2mm
Attachment	USB cable, user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP22, It means the device could protected against solid foreign objects of 12.5 mm and greater, and against vertically falling water drops when ENCLOSURE tilted up to 15°
Software version	V01
Device classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor charged Mode: Class II ME Equipment

WARNING: No modification of this equipment is allowed.

Complied European Standards List

Risk management	ISO/EN 14971:2012 Medical devices — Application of risk management to medical devices
Labeling	ISO/EN 15223-1:2012 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements
User manual	EN 1041: 2008 Medical equipment manufacturers to provide information
General Requirements for Safety	EN 60601-1: 2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance IEC/EN 60601-1-11: 2010 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment IEC/EN 80601-2-30:2009 Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated noninvasive sphygmomanometers
Electromagnetic compatibility	IEC/EN 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements	EN 1060-1:1995+A2:2009 Non-invasive blood pressure Part 1: General requirements EN 1060-3:1997+A2:2009 Non-invasive blood pressure Part 3: Supplementary requirements for electromechanical blood pressure measuring system
Clinical investigation	EN 1060-4: 2004 Automatic Blood Pressure Monitor overall system Interventional accuracy of the testing process
Usability	IEC/EN 60601-1-6: 2010 Medical electrical equipment -- Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability IEC/EN 62366: 2007 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	IEC/EN 62304:2006+AC: 2008 Medical device software - Software life cycle processes

EMC Guidance

1. The Blood Pressure Monitor needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS
2. Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance $d = 3,3$ m away from the equipment.

(Note. As indicated in Table 6 of IEC 60601-1-2:2007 for MEE EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3,3$ m at an IMMUNITY LEVEL of 3 V/m)

GUANGDONG TRANSTEK MEDICAL ELECTRONICS CO., LTD
Zone A, 5/F., Investment Building, No. 12, Huizhan East Rd., Torch
Development District, Zhongshan, Guangdong, 528437, China

MDSS - Medical Device Safety Service GmbH
Schiffgraben 41, 30175 Hannover, Germany

Terraillon SAS France & Headquarters

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Terraillon Asia Pacific Ltd

4/F, Eastern Centre
1065 King's Road
Quarry Bay - Hong Kong
Tel: + 852 (0)2960 7200
customerservice@terraillon-asia.com

Terraillon Corp USA

contact@terraillon.fr

Terraillon®

www.terraillon.com

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Terraillon UK Ltd