
# DAYVIA®

LUMINOTHERAPIE

## SlimStyle

FR	Notice d'utilisation	03
DE	Bedienungsanleitung	07
NL	Gebruiksaanwijzing	11
UK	User's guide	16

## INTRODUCTION

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Nederlands

## INTRODUCTION

Congratulations for purchasing this DAYVIA light therapy lamp.

Used in correct conditions, it will enable you to regain the light that you are lacking.

For certain people, the lack of light can cause various symptoms :

- irritability and sullen mood
- significant increase in appetite
- decrease of activity and the efficiency at work due to a lack of energy
- sleeping disturbances, chronic tiredness
- sexual disturbance (loss of libido)

Many international studies show that light therapy contributes to regaining energy, form, ability to concentrate and good mood, the same that you get whilst sunbathing.

The light therapy lamp also contributes to calming the sleeping rhythm for people who do night shifts, people who suffer from delayed sleep phase syndrome (falling asleep late and difficulties waking up in the morning) or advanced sleep phase syndrome (waking up too early) or people who travel and go across many time zones (jet-lag).

The effect of the light depends on its intensity but especially on its spectral curve (colour).

Indeed, some receptors located in our eyes are very sensitive to a light whose colour ranges between 460 and 480 nm (blue), they play a major part in our biological clock.

So a monochromatic blue light can prove to be as

efficient as a much more intense white light.

Therefore your device offers two possibilities which are as efficient as each other :

- an intense white light
- a monochromatic blue light

Please read this user's guide attentively in order to benefit from all the advantages of your DAYVIA lamp.

## SAFETY MEASURES

-Do not use the DAYVIA lamp in a damp environment (for example in a bathroom, next to a bathtub or shower cubicle, outdoors...).

-Keep the power cord away from heated surfaces.

-Only use the supplied DAYVIA power supply unit with your light therapy lamp. This power supply unit is specifically made for your lamp, it must not be used with any other equipment.

-The plug serves as the power switch. The device must be used with an accessible plug.

-This lamp must not be used by children, or physically or mentally disabled people.

-Any modifications to this lamp are forbidden.

-The light must not be used if the cord and/or the power supply unit are apparently damaged or defective.

## PRECAUTIONS

The light therapy lamp must be used in an environment lit by another source of light, the lamp is complementary to normal light (electric or natural light).

Contrary to natural daylight or to other lamps on the market, this light therapy medical device diffuses no UV rays.

Studies have shown that light therapy can be used in safety (for example your light therapy lamp diffuses a specific blue light or an intense white light – up to 10 000 lux for the white one, whereas on a beach in the summer, the light can reach up to 100 000 lux as well as harmful rays such as UV of which your lamp is devoid).

Nevertheless, before starting your treatment, please consult your general practitioner if :

- You suffer or have suffered from depression
- You are hypersensitive to light
- You suffer from an eye illness
- You suffer from a mental condition
- You are undergoing a treatment with psychotropic or photo-sensitizing medication

If you encounter any problems (for example headaches or eye tiredness) after the use of the light therapy lamp, please tell your general practitioner. Often, all you have to do is start the treatment with shorter sessions (10 minutes) and gradually increase them each day until the recommended duration.

## TURNING ON/OFF

The Dayvia Slim Style lamp is carefully controlled before packaging and expedition.

You will find the user's guide, the power supply unit and your lamp in the package.

We recommend you keep the packaging, in case return is necessary, else it must be disposed of in accordance with your country's regulations.

Here are the different steps to turn on your DAYVIA light therapy lamp :

1. Insert the male connector of the lamp onto the female connector situated on the supplied power supply unit.
2. Check that the voltage indicated on the power supply unit corresponds to the local sector voltage, then plug in the device.
3. Turn on the lamp, in a lit environment, using the on/off  switch located at the top of the lamp (when the lamp is in the vertical position).

You can then adjust the intensity with the + and – buttons next to the on/off switch.

When you press on the switch once, you will get an intense white light which is very efficient from a therapeutic perspective. You can reduce its intensity with the – button, the light will become warmer and warmer for back-up lighting. When you press on the on/off switch a second time, you will get a low-intensity monochromatic blue light which is as efficient as the

English

English

intense white light.

- When finished using the device, turn it off by simply pressing the on/off switch.
- Unplug the device from the socket.

## USER'S RECOMMENDATIONS

Place the lamp between 30 and 70cm from your eyes in the field of vision (at an angle of approximately 30° in front of your face), it is not necessary to continually look at the lamp, pursue your activities normally (breakfast, office work, computing, make-up...).

For a pleasant use, without a dazzling effect, make sure that the room is lit (especially for the first days of use) to minimize the contrast between the lighting of the lamp and that of the environment.

It is advised to start the light therapy sessions as soon as you feel a lack of light, before the first symptoms of seasonal affective disorder (SAD) appear.

As a prevention or treatment for SAD, it is advised to carry out sessions in the morning after waking up.

It is also advised to do one session every morning in 3-week cycles during the dark periods of the year. These figures are to be adapted to each person depending on the results after a few days of treatment.

### Maximum intensity white light option:

### Maximum intensity blue light option:

Generally, all people wanting to compensate their lack of light (SAD, working in a confined room...) with light therapy must give priority to morning exposure : we recommend an exposure of 30 minutes to 1 hour (see above for the distance) at maximum intensity. Morning exposure, ideally just after waking up is recommended for people with delayed sleep phase syndrome (falling asleep late and difficulties waking up in the morning) and for people preparing a trip easterly (crossing more than 3 time zones), the effects of jet-lag will be notably minimized.

Evening light therapy is reserved for people working at night, people with advanced sleep phase syndrome (waking up too early) and for people preparing a trip westerly (crossing more than 3 time zones) to avoid jet-lag.

After your daily light therapy session, you can use your Dayvia Slim Style lamp all day long as a back-up light or a desk lamp by reducing the intensity of the white light according to your wishes.

## DEVICE MAINTENANCE

Unplug the lamp and the power supply unit and leave it to cool down before cleaning it.

Never use abrasive substances or detergents to clean the device : use a damp cloth making sure that no liquid penetrates into the device.

Being a LED technology, there is no maintenance and the lifespan is over 30 000 hours.

English

## TECHNICAL CHARACTERISTICS

The characteristics are prone to modifications without notice.

- Power supply : 21W
- Power supply input voltage (V AC) : 220-240 V
- Power supply output voltage (V DC) : 6 V
- Frequency : 50-60 Hz
- Insulation class : II
- Altitude of use < 2000m
- Temperature of use : from +10°C to +40°C
- Relative humidity level whilst in use : from 30% to 75%
- Storage temperature : from -5°C to +45°C
- Over voltage category II
- Length of cords : 3.2 m
- Dimensions of the device : 34.6 x 13.7 x 9.6 (h x w x d)
- Device weight : : 1.1 kg
- Device weight with power supply : 1.3 kg
- Class II A medical device : CE0459

At the end of numerous trials in 3 French independent laboratories accredited COFRAC, DAYVIA slim style is certified conform to all the guideline norms 93/42/CEE and 2007/47/CEE applied to medical devices in Europe :

- conforms to electrical security
- conforms to EMC (electro magnetic fields)
- conforms to photobiological risks of light : group 0 (no risk).

## THE SYMBOLS ON THE POWER SUPPLY UNIT

Consult the supplied user's guide.

This device is class II.

This device must not be discarded with household wastes, take it to a waste disposal facility where it can be recycled. By doing so, you will contribute to the protection of the environment.

Equipment with safety insulation.

<sup>0459</sup> This device is certified conform to the guideline regulating medical devices.

**DAYVIA®**  
LUMINOTHERAPIE

Immeuble Le leeds  
253 Bd du Leeds  
59777 LILLE

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0459 D-Notice Slim Style - Rev B - septembre 2014